



Happier Living Everyday

BPL Medical
Technologies



**Together in your fight
against Covid 19**

Follow Us on





Critical Care

Patient Monitors



BPL Excello

- High Resolution Display
- Standard Configuration Standard configuration
- Calculation
- Arrhythmia and ST Analysis
- 168 Hour Trend review
- Waveform Display
- Connectivity to CNS
- Optional Parameters & Configuration



BPL Vivid Vue

- High Resolution Display
- Enhanced Connectivity Solution
- Standard Configuration
- High performance Li-ion Battery
- Optional Parameters & Configuration
- Comprehensive Data Storage

Oxygen Concentrators



BPL Oxy 5 Neo | Oxy 5 Neo Dual

- **Oxy 5 Neo:** Single outlet | **Oxy 5 Neo Dual:** Double outlet
- Flow range: 0.5 to 5L/min; 95.5% O₂ concentration
- Built-in nebulizer
- LCD Screen to view usage hours
- LED Lights for System Status, Audio alarm for power failure, Low Oxygen concentration, compressor failure, pressure cycle failure
- Timer Function for auto shut down
- Flame retardant thermo plastic molded cabinet for strength and durability
- Storage compartment to place accessories
- Portable, user-friendly, safe, quiet and reliable

Video Laryngoscope



BPL Video Laryngoscope

- High Resolution 3.5" Full View LCD Monitor
- 66° Horizontal Screen Rotation along with Vertical movements for enhanced clinical assistance
- 2MP CMOS camera, High Definition Anti Fog lens and integrated LED light source for high performance visualisation
- Videos and images can be transferred via a micro USB port for easy storage and access
- Rechargeable Lithium ion battery with a 200 minutes capacity with a low battery indicator
- 4GB built in memory, for continuous video recording and image capturing with an use of a button.
- Reusable Blades

ECG

Powered by
Glasgow ECG Interpretation
Algorithm



Cardiart 9108D

- Twelve Channel ECG Recording
- 7-inch High resolution color LCD with touch screen option
- Glasgow ECG Interpretation Algorithm
- Internal storage 800 ECGs and can be enlarged by USB & SD card
- Li-ion rechargeable battery – 300 ECG's
- Report formats: PDF and optional SCP/FDA-XML/DICOM data export
- Connectivity: FTP, Supports barcode scanner & USB Printer, optional Wi-Fi
- Optional PC interface software allows to access EMR/ PACs/ GDT-based Systems (HL7 protocol)

Powered by
Glasgow ECG Interpretation
Algorithm



Cardiart GenX3

- Three Channel ECG Recording with trace darkness control
- 4.3-inch color LCD screen
- Glasgow ECG Interpretation Algorithm
- Alphanumeric keypad with soft silicone keys
- Full disclosure view of ECG analysis on display
- Internal storage 250 ECGs and can be enlarged by USB
- Li-ion rechargeable battery – 250 ECG's
- Connectivity: supports USB Printer, optional FTP, Export selected record(s) in HL7 format
- Optional: Real-time ECG transfer to PC through USB

Syringe & Volumetric Pump



Acura S1

- Dual Display: 2.8-inch color LCD & LED for rate display
- Compatible with all standard syringes: 5ml – 60ml
- Infusion Modes: Rate Mode, Time Mode and Drug/Weight Mode
- In-built Drug Library
- Dual Bolus Modes: Manual & Automatic bolus
- 3 level occlusion sensitivities setting (programmable)
- 2000 history records and last configuration memory
- Horizontal & Vertical fixation clamp



Acura V1

- 3.2-inch High definition color display
- Infusion Modes: Rate Mode, Time Mode, Weight Mode, Stage Mode, Drip Rate Mode
- In-built Drug Library
- Dual Bolus Modes: Manual & Automatic bolus
- 3 level of occlusion pressure (programmable)
- Compatible with all brands of infusion sets with IV set library
- 2000 history records and last configuration memory
- Horizontal & Vertical fixation clamp



ICU Ventilators



Prisma Vent40

- Quiet operation with peak flow of 220 Litres.
- Caters ET Tube and Mask application to Pediatric and Adult Pa-tients.
- Pressure control and spontaneous modes to meet patient de-mand.
- Tidal volume from 100ml – 2000ml can be targeted to assure vol-ume delivery
- Volume & Flow waveforms can be viewed & Alarms can be man-aged
- Low flow oxygen up to 30 Liters can be connected with machine.



Elisa 800 VIT/600

- **ELISA 800VIT/600:** 18.5" large display helps in better monitoring of patients
- Trolley mounted premium class of ventilators suitable for Adult, Paediatric and Neonatal applications with highest quality of ventilation
- Capable of ventilating in both Invasive and non-invasive modes for adults and children with a weight of more than 2 Kgs
- Instant View Technology that gives the operator an intuitive grasp of the patient's situation
- Optional Upgradation to automatic cuff pressure management and transpulmonary pressure measurement & VIT (Ventilator Integrated Tomography)



Imaging

Mobile X-ray Solutions



BPL MRAD[®] 3.5 DR

- Powered by Imported X-ray generator for proven quality and reliability
- 3.5kW, 110kV Output with a 40KHz high frequency X-ray generator
- 1.5mm² Focal spot produces sharper images
- User friendly X-Ray control panel with alphanumeric LCD display & 155 APR
- Concealed cable design in Arm for hygienic advantages in hospital use
- Wireless 14"X17" Flat Panel detector with Directly deposited CsI
- Detector with IPX-4 rating for ingress protection from liquids
- 12.3' Touch console for software provides clinicians access to high quality images
- PACS connectivity using LAN/Wi-Fi aids in easy access and transfer of images

Approved By  



BPL MRAD[®] 3.6+ DR

- 110kHz ultra high frequency technology ensures that only useful dose is produced over the entire duration of the exposure while suppressing the skin dose.
- 3.5kW, 110kV output X-Ray generator.
- 1.5mm² Focal spot produces sharper images
- User friendly X-Ray control panel with alphanumeric LCD display & 155 APR
- Concealed cable design in Arm for hygienic advantages in hospital use
- Wireless 14"X17" Flat Panel detector with Directly deposited CsI
- Detector with IPX-4 rating for ingress protection from liquids
- 12.3' Touch console for software provides clinicians access to high quality Images
- PACS connectivity using LAN/Wi-Fi aids in easy access and transfer of images

Approved By  



BPL MRAD[®] 3.6+


- 110kHz ultra high frequency technology ensures that only useful dose is produced over the entire duration of the exposure while suppressing the skin dose.
- 3.5kW, 110kV output X-Ray generator.
- Stationary anode X-Ray tube with 1.5mm² smaller focal spot produces sharper images
- Concealed cable design in Arm for hygienic advantages in hospital use.
- 155 anatomical programming feature for precise selection of KV & mAs based on patient's anatomy (Optional)
- Ergonomically design product with easy maneuverability & more storage space

Approved By  



BPL MRAD[®] 100

- 100KV, 100mA Line frequency mobile X-Ray
- Alphanumeric LCD display with 4 X 20 characters for easy viewing
- Precise selection of KV (2KVp / Step)
- 155 anatomical programming feature for precise selection of KV & mAs based on patient's anatomy
- Ergonomically design product with easy maneuverability & more storage space



Ultrasound Solutions



E-CUBE 8 DIAMOND[#]

- Latest i7 processor with SSD for high frame rate images.
- Transducers powered with PowerView™ to dissipate heat for stable images.
- A 21.5" high-definition LED main monitor and a 10.4" LED touch screen defined for convenience.
- CUBE Strain™ and TDI available for advance cardiac applications
- Advanced user interface with power presets to reduce key strokes.
- High performance linear transducers and 230° FOV for endovaginal images.
- Articulating Arm
- Control panel height adjustment



E-CUBE i7[#]

- A portable machine built around a cart based system architecture with 15.6" LED screen.
- Optimal imaging suite with needle vision.
- Echocardiography available with TDI.
- Crystal signature™ (a single crystal technology) and high frequency probes available.
- A sliding keyboard cover, the first of its kind, reducing the stress on wrists.
- Tripple Connector Trolley



Home Care



Pulse Oximeters



Nebulizers



**Non-Contact
Thermometers**



BP Monitors



CPAP/ BiPAP



AIR PURIFIERS



PPE



Protective Masks



BPL Goggles



BPL Face shield



BPL PPE Kit



Coverall with hood



Head Cover



3Ply mask



Latex Gloves



Shoe cover



**Always closer
to You**

HEAD OFFICE

BPL Medical Technologies Private Limited
11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076
Ph: +91 80 26484388/ 2648 4350
Toll Free: 1800-4252355
Website: www.bplmedicaltechnologies.com
For Enquiries: sales.medical@bpl.in

1800-425-2355


9 AM to 8 PM (Monday-Saturday)

SALES OFFICES

Ahmedabad	079-40031471
Bangalore	080-26484350/ 26484388
Bhopal	0755-4059311
Chandigarh	0172 - 2970883
Chennai	044-24353074/ 24353013
Cochin	0484-2558418 / 4014406
Hyderabad	040 - 27017175

Jaipur	0141 - 2525933
Kolkata	033 - 24434071/ 24434072
Lucknow	0522-4042194/ 4065350
Mumbai	022 - 61046000
Delhi	011-29248313/9248314/29248315
Factory	0491 - 2573100 / 2573105

A Global Medical Technology Company



Service Helpline
1800-425-2355

CERTIFIED ISO 13485 : 2016 COMPANY

BPL MEDICAL TECHNOLOGIES PRIVATE LIMITED

Regd. Office: 11th KM, Bannerghatta Road,
Arakere, Bangalore - 560076, India.

Ph: +91 80 26484388/ 2648 4350

Toll Free: 1800-4252355

Website: www.bplmedicaltechnologies.com

For Enquiries: sales.medical@bpl.in

CIN: U33110KA2012PTC067282



Happier Living Everyday