

# EXCELLO



CE 0843

MULTI PARAMETER MONITOR



*Happier Living Everyday*

# Experience the Finest Display



# Experience the Finest Display

CE 0843



 The device provides arrhythmia detection, ST & HRV analysis with 3-lead/5-lead ECG monitoring and comprehensive trend data review



#### High resolution Display

High resolution 10.4" colour TFT display



#### Standard Configuration

Standard configuration of ECG, RESP, SpO<sub>2</sub>, NiBP, TEMP



#### Drug Dose Calculation

The added functions of the device are drug dose calculation /HRV analysis/OXY CRG



#### Arrhythmia and ST Analysis

The device provides for arrhythmia detection and ST analysis



#### ECG Monitoring

3-Lead/5-Lead monitoring



#### Waveform Display

Maximum 10 waveforms display



#### 168 Hour Trend

Provides for maximum of 168 hours of trend review



#### Pacemaker Detection

Detects pacemaker pulses



#### Rechargeable Battery

>120 minutes operating time in fully charged condition and >240 minutes with additional battery (optional)



#### Networking

Wired or wireless networking with Central monitoring system software (optional)



#### Optional Parameters & Configurations

12 Lead ECG, 2-TEMP, 2-IBP, ETCO<sub>2</sub> (mainstream/ side stream/ micro stream), AGM, ICG, thermal recorder and lithium battery