

DIGITAL BLOOD PRESSURE MONITOR

BPL
home care

A fully automatic digital blood pressure monitor with 2.5" LCD display. Works on oscillometric principle and automatic pressurization technology reduces measurement time and provides accurate readings. Dual user memory function, each result will be displayed on the screen and stored automatically. WHO BP classification indicator helps user for easy detection of blood pressure level. Device comes with latex free adult cuff which is safe and convenient for user. Also works on Micro USB Cable charger.

BPL 120/80 B16



SPECIFICATIONS

- **DISPLAY**- 2.5" LCD DISPLAY
- **MEASUREMENT RANGE**- NIBP: 0 - 299mmHg PULSE RATE: 40 - 180BPM
- **ACCURACY**- ± 3 mmHg / $\pm 5\%$ of Reading
- **CUFF SIZE**- 22 - 32CM, RECOMMENDED FOR ADULTS
- **INTERPRETATION**- WHO COLOUR CODED BP CLASSIFICATION INDICATOR
- **MEMORY**- AUTOMATIC MEMORY FUNCTION RECORDS THE LAST 99 RECORDINGS IN 2 SEPARATE USER ACCOUNTS FOR EASY TRACKING
- **IHB DETECTION**: ALERTS USER IN CASE OF IRREGULAR HEARTBEAT DETECTION
- **AUTO SHUT OFF** - AUTOMATIC SHUTDOWN FUNCTION ENHANCE BATTERY LIFE
- **LOW BATTERY INDICATION** - VISUAL INDICATION ON LCD FOR LOW BATTERY
- **POWER SOURCE** - 4 x AA BATTERY OR DC5V, 1A MICRO USB ADAPTER

BPL ORDERING CODE: 91MED300

BOX CONTENTS- Unit, Cuff 22 - 32cm, 4 x AA batteries, Micro USB cable & User manual

CONSUMABLES- Cuff 22 - 32cm, 4 X AA battery, Miro USB cable

BPL Medical
Technologies